


EnligHTN™

RF Ablation Generator



AN ENLIGHTENED APPROACH TO MULTI-ELECTRODE RENAL DENERVATION

The EnligHTN radiofrequency (RF) ablation generator is the power source for the EnligHTN renal denervation system. The generator controls the delivery of consistent transmural ablations, enabling effective renal denervation.^{1,2} The temperature-controlled system and default pre-set settings make it easy to use, with features including temperature control, pre-ablation diagnostics and automatic shutoff.

PROPRIETARY, TEMPERATURE-CONTROLLED ALGORITHM

- Pre-programmed settings for use with the EnligHTN multi-electrode renal artery ablation catheter.
- Uses a proprietary, temperature-controlled algorithm developed to produce transmural lesions for effective denervation.^{1,2}

INDEPENDENT CONTROL OF ELECTRODES

- Each of four monitoring channels can be independently monitored and adjusted for enhanced control of the ablation process.
- Electrodes can also be selectively disabled if desired.

SAFETY SHUTOFF

- Generator automatically shuts down at the completion of ablation, or in the event temperatures or impedance levels exceed pre-set limits.



ST. JUDE MEDICAL™
MORE CONTROL. LESS RISK.

YOUR PARTNER IN RENAL DENERVATION TREATMENT

The EnligHTN renal denervation system offers advanced renal artery ablation treatment and is designed to deliver radiofrequency (RF) energy to the renal nerves to achieve targeted denervation. The core system components are the EnligHTN RF ablation generator, the EnligHTN renal artery ablation catheter and the EnligHTN guiding catheter.

Ordering Information

EnligHTN Renal Denervation System	
Reorder Number	Product Name
ENL-T115	EnligHTN RF Ablation Generator
ENL-8F55	EnligHTN Guiding Catheter
ENL-LG-01	EnligHTN Renal Artery Ablation Catheter - large
ENL-SM-01	EnligHTN Renal Artery Ablation Catheter - small
ENL-CC-01	EnligHTN 1641 Connector Cable
PREP-01	Patient Return Electrode
ENL-PC-01	EnligHTN Power Cable - replacement



EnligHTN RF Ablation Generator



EnligHTN Guiding Catheter



EnligHTN Renal Artery Ablation Catheter

1. St. Jude Medical. Data on file. 900900781, 2010-2011 Pre Clinical Studies Summary.
2. St. Jude Medical Instruction for Use, Ablation catheter #10006654702, Generator #10006927602.

ATRIAL FIBRILLATION

Global Headquarters
One St. Jude Medical Drive
St. Paul, Minnesota 55117
USA
+1 651 756 2000
+1 651 756 3301 Fax

St. Jude Medical Brasil Ltda.
Rua Itapeva, 538
5º ao 8º andar
01332-000 – São Paulo – SP
Brazil
+55 11 5080 5400
+55 11 5080 5423 Fax

CARDIAC RHYTHM MANAGEMENT

Cardiovascular Division
177 East County Road B
St. Paul, Minnesota 55117
USA
+1 651 756 4470
+1 651 756 4466 Fax

St. Jude Medical Japan Co., Ltd.
Shiodome City Center 15F
1-5-2, Higashi-Shinbashi
Minato-ku
Tokyo 105-7115
Japan
+81 3 6255 6370
+81 3 6255 6371 Fax

CARDIOVASCULAR

SJM Coordination Center BVBA
The Corporate Village
Da Vincilaan 11 Box F1
1935 Zaventem, Belgium
+32 2 774 68 11
+32 2 772 83 84 Fax

NEUROMODULATION

St. Jude Medical (Hong Kong) Ltd.
Suite 1608, 16/F Exchange Tower
33 Wang Chiu Road
Kowloon Bay, Kowloon
Hong Kong SAR
+852 2996 7688
+852 2956 0622 Fax

SJMenligHTN.com
SJMprofessional.com



ST. JUDE MEDICAL™
MORE CONTROL. LESS RISK.

Product referenced is approved for CE Mark.

EnligHTN, ST. JUDE MEDICAL, the nine-squares symbol and MORE CONTROL. LESS RISK. are registered and unregistered trademarks and service marks of St. Jude Medical, Inc. and its related companies. ©2012 St. Jude Medical, Inc. All rights reserved.

Item GMVASC401EN IPN 2156-12